


# In vitro Proof of Principle Experiment with the Osmotic Pressure-based Sencell Implantable Glucose Sensor Technology

Ramljak S., Lecher B., Meier A., Van Wyk E., Watson S., Frisvold R., Pfützner A., Pfützner Science & Health Institute and Sciema, Mainz, Germany; MfD, Wendelsheim, Germany; Cambridge Consultants, Cambridge, UK; Lifecare Bergen, Norway

## Background

Minimally invasive or implantable continuous glucose monitoring systems have become available and are increasingly adopted. A novel implantable glucose affinity biosensor (Sencell, LifeCare, Norway) is based on competitive and reversible binding of glucose and polysaccharide dextrane to the glucose specific lectin concavalin A (ConA). The Sencell technology uses osmotic pressure difference arising between a reagent chamber (containing active fluid with ConA and dextrane) and a diffusion chamber (in direct contact with interstitial fluid) to determine interstitial glucose concentrations. Both chambers are separated by a nanoporous membrane permeable to glucose and water, but not to ConA or dextrane. Because the skin anatomy of pigs markedly resemble the human situation, we performed a first proof-of-principle experiment in this animal model using a wired device prototype.

Fig.1.: Osmotic pressure cell principle and experimental setup





## Methods

Implantation of four Sencell sensors and one Dexcom® G4 control sensor was performed in the back and the neck area of three pigs (one female, two male) respectively. Subsequently, they received an oral glucose load of 75 to 150 g of dextrose. Reference measurements from capillary blood samples were performed using the YSI 2300 STAT Plus glucose analyzer every 15 min for 5 h. Statistical analysis tried to identify sensor patterns from the crude pilot sensor prototypes indicating the ability of the osmotic sensing technology to track interstitial glucose

## Results

Several Sencell prototype devices were able to track glucose changes, especially when blood glucose levels exceeded 200 mg/dL. The magnitude of the signals was in a predicted range and working sensors were matching G4 results and displayed a lag phase of 20 minutes vs. plasma reference glucose. Although the body temperatures in the animals within three experimental days partially exceeded 39°C, the activity of ConA was preserved. No external signs of inflammation were observed in the histology examinations.

Fig.2.: Raw and retrospectively calibrated Sencell signals in comparison to the different reference methods



## Conclusions

In conclusion, Sencell proof-of-principle has been demonstrated *in vivo*. Next development steps are now targeting to obtain a working prototype with more stable and predictable performance.