

VeraPrep Biotin™

Part No. 400014

For Research Use Only

NAME OF THE PRODUCT

VeraPrep Biotin

INTENDED USE

VeraPrep Biotin is a research use only sample pre-treatment reagent that removes biotin interference from a serum or plasma sample. It is used in conjunction with a non-treated sample to investigate possible interference from biotin in immunoassays.

VeraPrep Biotin is not intended for use in diagnosis of disease or other conditions, including a determination of the state of health, in order to cure, mitigate, treat, or prevent disease or its sequelae.

SUMMARY AND EXPLANATION

Biotin, also known as vitamin B7, is a water-soluble B vitamin often found in multi-vitamins and over the counter health and beauty supplements. *In vitro* laboratory tests that employ streptavidin-biotin binding mechanisms have the potential to be affected by high circulating biotin concentrations. Biotin can be attached through covalent bond to a variety of targets—from large antibodies to steroid hormones—with minimal effect on their specific non-covalent binding with avidin, streptavidin, or NeutrAvidin proteins. Therefore, biotin has been frequently used in the detection systems of immunoassays of different forms.

Immunoassays are generally categorized as either sandwich immunoassays (non-competitive) or competitive inhibition immunoassays. In general, streptavidin-biotin binding is used during assay incubation to couple biotinylated antibodies in sandwich immunoassays, or biotinylated antigens in competitive immunoassays, to streptavidin-coated surfaces. When a biological specimen contains excess biotin, the biotin competes with the biotinylated antibodies or antigens for binding to the streptavidin-coated surfaces, resulting in reduced capture of the biotinylated antibodies or antigens. Excess biotin produces falsely low results in sandwich immunoassays because the assay signal is directly proportional to the analyte concentration. Excess biotin in competitive immunoassays causes falsely elevated results because the assay signal is inversely proportional to the analyte concentration. Specific details of biotin interference have been extensively described in other publications.(1-13)

Normal circulating concentrations of biotin derived from the diet and normal metabolism are too low (< 1 ng/mL) to interfere with biotinylated immunoassays. However, ingestion of high-dose biotin supplements (e.g., 5 mg or higher) can result in significantly elevated blood concentrations that can interfere with commonly used biotinylated immunoassays. In certain medical conditions, extremely high biotin doses (e.g., 100 mg or higher) can result in serum or plasma biotin levels of > 1000 ng/mL.(14)

According to the FDA, biotin in blood or other samples taken from patients who are ingesting high levels of biotin



can cause falsely high or falsely low results in biotin-based immunoassays, depending on the design of the assay.(14)

Biotin interference thresholds differ widely among assays, even on a single platform. Tests with biotin interference thresholds < 51 ng/mL are considered high risk tests, or vulnerable immunometric and competitive methods.(1)

VeraPrep Biotin is a research use only sample pre-treatment reagent that can be used to help rule-in or rule-out biotin interference. It uses a 15 minute procedure to remove up to 500 ng/mL free biotin in serum or plasma without sample dilution, and it can be used to remove biotin levels > 500 ng/mL using an enhanced procedure.

REAGENT

VeraPrep Biotin reagent contains proprietary superparamagnetic nanoparticles covalently conjugated to Streptavidin. After > 30 seconds magnetic separation using VeraMag™ (Part No. 40020), the VeraPrep Biotin storage buffer is aspirated and discarded, the serum or plasma sample is added and mixed, and the reagent incubates with the sample to bind and remove biotin interference from the sample. After a 10 minute incubation, the reagent is magnetically separated for > 4 minutes using VeraMag and the sample supernatant is aspirated and saved for testing. Each VeraPrep Biotin vial contains enough reagent to pre-treat 20 different 400 µL samples using the standard procedure. The reagent is in the form of a liquid and must be well mixed prior to use to ensure homogeneous resuspension of the nanoparticles.

REAGENTS AND MATERIALS PROVIDED

Streptavidin coated superparamagnetic nanoparticles in TRIS buffer and detergent. Preservative: 0.05% sodium azide.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Pipetting device(s) capable of delivering 50 µL up to 1000 µL
2. Disposable pipette tips
3. Micro tube 2ml with cap (SARSTEDT Order Number 72.694)
4. Vortex mixer
5. VeraMag (Part No. 40020)
6. Timer
7. Laboratory mixer
8. Transfer tube
9. Personal protective equipment

STORAGE AND STABILITY

Upon receipt, store in the box at 2° - 8°C. The shelf life is approximately six (6) months. Refer to the expiration dates marked on the vial label.

FOR RESEARCH USE ONLY

WARNINGS AND PRECAUTIONS

1. For Research Use Only. Not for use in diagnostic procedures.
2. Do not use test components beyond their expiration dates.
3. This product contains sodium azide. For a specific listing, refer to the **REAGENTS AND MATERIALS PROVIDED** section. This material and its container must be disposed of in a safe way.
4. Dispose of all potentially contaminated test components in a biohazard container.
5. Each box contains 1 foam vial holder (donut) to hold the VeraPrep Biotin vial during use and to prevent it from accidentally falling over and spilling reagent.
6. Remove the reagent storage solution using VeraMag before adding the sample to prevent sample dilution.
7. VeraPrep Biotin should be used with SARSTEDT tubes (Order Number 72.694). Other tubes types have not been studied.
8. Do not incubate the VeraBind Biotin reagent on VeraMag without any storage solution or sample.

SPECIMENS COLLECTION AND PREPARATION

Follow manufactures specification for blood collection and serum or plasma preparation.

STANDARD PROCEDURE:

VeraPrep Biotin has a binding capacity ≥ 1 ng biotin per 1 μ L reagent.

The standard VeraPrep Biotin procedure uses a 1:2 ratio of VeraPrep Biotin reagent to serum or plasma sample, or 200 μ L reagent and 400 μ L sample, to deplete biotin interference up to 500 ng/mL. Smaller and larger sample volumes can be used if a 1:2 ratio of reagent:sample is maintained.

Standard Procedure 1:2 Ratio (Reagent:Sample)		
VERAPREP Biotin (μ L)	Serum or Plasma (μ L)	Samples (Uses per Vial)
50	100	80
100	200	40
200	400	20
400	800	10

VeraPrep Biotin sample pre-treatment instructions to remove biotin interference ≤ 500 ng/mL from serum or plasma:

1. Remove the VeraPrep Biotin reagent vial from storage and vortex for a minimum of 10 seconds at medium speed to mix well and resuspend the reagent.
2. Insert the reagent vial in the foam vial holder.
3. Insert an empty Micro tube 2ml (SARSTEDT Order Number 72.694) into the VeraMag magnet until the collar of the tube contacts the magnet frame.
4. Dispense **200 μ L** of the well-mixed **reagent** into the empty tube to separate the reagent on the magnet for > 30 seconds to form a reagent pellet.
5. Carefully aspirate and discard all of the storage buffer supernatant (~ 200 μ L) without disturbing the reagent pellet.

6. Dispense **400 μ L** of well-mixed serum or plasma **sample** into the tube containing the reagent pellet.
7. Tighten the screw cap on the tube, remove the tube from the magnet, and vortex for a minimum of 10 seconds at medium speed to mix well and resuspend the reagent in the sample.
8. Place the tube onto a laboratory mixer at medium speed and incubate at room temperature for 10 minutes.
9. Loosen and remove the screw cap and insert the tube into the magnet until the collar of the tube contacts the magnet frame.
10. Magnetically separate the reagent for > 4 minutes to form a reagent pellet.
11. Carefully aspirate the sample supernatant without disturbing the reagent pellet and dispense the sample into a transfer tube for testing. Note: All of the sample supernatant (~ 400 μ L) can be aspirated if this step is performed carefully. If any of the reagent is accidentally aspirated then simply return the sample/reagent mixture to the tube and return to step 10.
12. The sample is now ready for testing.

ENHANCED PROCEDURE

Extremely high biotin doses (e.g., 100 mg or higher) can result in serum or plasma biotin levels of > 1000 ng/mL. (14)

The enhanced VeraPrep Biotin procedure uses a 3:2 ratio of VeraPrep Biotin reagent to serum or plasma sample, or 600 μ L reagent and 400 μ L sample, to deplete biotin interference greater than 500 ng/mL. Smaller and larger sample volumes can be used if a 3:2 ratio of reagent:sample is maintained.

Enhanced Procedure 3:2 Ratio (Reagent:Sample)		
VERAPREP Biotin (μ L)	Serum or Plasma (μ L)	Samples (Uses per Vial)
150	100	26
300	200	13
600	400	6
1200	800	3

VeraPrep Biotin sample pre-treatment instructions to remove biotin interference > 500 ng/mL from serum or plasma:

1. Remove the VeraPrep Biotin reagent vial from storage and vortex for a minimum of 10 seconds at medium speed to mix well and resuspend the reagent.
2. Insert the reagent vial in the foam vial holder.
3. Insert an empty Micro tube 2ml (SARSTEDT Order Number 72.694) into the VeraMag magnet until the collar of the tube contacts the magnet frame.
4. Dispense **600 μ L** of the well-mixed **reagent** into the empty tube to separate the reagent on the magnet for > 30 seconds to form a reagent pellet.
5. Carefully aspirate and discard all of the storage buffer supernatant (~ 600 μ L) without disturbing the reagent pellet.
6. Dispense **400 μ L** of well-mixed serum or plasma **sample** into the tube containing the reagent pellet.

- Tighten the screw the cap on the tube, remove the tube from the magnet, and vortex for a minimum of 10 seconds at medium speed to mix well and resuspend the reagent in the sample.
- Place the tube onto a laboratory mixer at medium speed and incubate at room temperature for 10 minutes.
- Loosen the screw cap and insert the tube into the magnet until the collar of the tube contacts the magnet frame.
- Magnetically separate the reagent for > 4 minutes to form a reagent pellet.
- Carefully aspirate the sample supernatant without disturbing the reagent pellet and dispense the sample into a transfer tube for testing. Note: All of the sample supernatant (~ 400 µL) can be aspirated if this step is performed carefully. If any of the reagent is accidentally aspirated then simply return the sample/reagent mixture to the tube and return to step 10.
- The sample is now ready for testing.

LIMITATIONS OF USE


- For Research Use Only. Not for use in diagnostic procedures.
- VeraPrep Biotin is not intended to replace manufacturer controls provided with the primary assay.
- The sample may need to be submitted for further research for potential false negative or false positive assay results due to human anti-Streptavidin interference. (15-20)

PERFORMANCE CHARACTERISTICS


An internal study was completed to demonstrate the ability of VeraPrep Biotin to remove high levels of biotin from serum samples and to rule-in biotin interference in a PTH ELISA test susceptible to biotin interference (i.e. Streptavidin coated microtiter plate wells, sample added to microtiter plate wells prior to the biotinylated capture antibody).

- Four apparently healthy adult volunteers had fasting baseline serum samples collected, ingested 100-200 mg of over the counter (OTC) biotin, and had serum samples collected again 1 hour post-biotin ingestion. All serum samples were aliquoted and stored frozen at -80° Celsius until testing.
- The baseline serum samples were tested by the PTH Intact ELISA (DRG PTH Intact ELISA, Part No. EIA-3645). PTH values ranged from 28.1 to 50.3 pg/mL.
- The 1 hour post-biotin ingestion samples were tested by the PTH Intact ELISA. All PTH results were < 1.57 pg/mL, or below the lower limit of detection (LLD).
- Endogenous biotin levels in each sample were determined by LC-MS/MS (Department of Lab Medicine, University of Washington Medical Center, Seattle, WA). Biotin concentrations ranged from 294 to 861 ng/mL.

- Samples 1 and 4 had biotin levels < 500 ng/mL and were pre-treated using the VeraPrep Biotin standard procedure.
- Samples 2 and 3 had biotin levels greater than 500 ng/mL and were pre-treated using the VeraPrep Biotin enhanced procedure.
- To verify biotin interference removal, VeraPrep Biotin pre-treated samples were tested using the Immundiagnostik IDK® Biotin ELISA kit (Part No. K8141, measuring range of 48.1 – 1,100 ng/L). Biotin concentrations ranged from 0.2 to 1.0 ng/mL, or within normal plasma levels (200 – 1,200 ng/L).



- Immediately after VeraPrep Biotin pre-treatment of Samples 1-4, PTH values were measured using the PTH Intact ELISA. PTH values ranged from 26.7 to 52.0 pg/mL.
- The 1 hour post-biotin ingestion samples had high levels of biotin interference per LC-MS/MS (294 to 861 ng/mL) and undetectable PTH values by the PTH Intact ELISA (< 1.57 pg/mL).
- The 1 hour post-biotin ingestion samples pre-treated with VeraPrep Biotin had physiologically normal biotin values per the Biotin ELISA kit (< 1.1 ng/mL) and normal PTH values by the PTH Intact ELISA (26.7 to 52.0 pg/mL).



- When comparing VeraPrep Biotin treated samples (grey bars) to the baseline PTH values (blue bars) the results recovered from 95 to 113% (mean recovery of 105%).

REFERENCES

1. Samarasinghe S, Meah F, Singh V, Basit A, Emanuele N, Emanuele MA, Mazhar Ai, Holmes EW. Biotin interference with routine clinical immunoassays: understand the causes and mitigate the risks. *Endocrine Practice*: August 2017, Vol. 23, No. 8, pp. 989-998.
2. Grimsey P, Frey N, Bendig G, Zitzler J, Lorenz O, Kasapic D, Zaugg CE. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *International Journal of Pharmacokinetics* 37. Online publication date: 14-Sep-2017.
3. Samarasinghe S, Meah F, Singh V, Basit A, Emanuele N, Emanuele MA, et al. Biotin interference with routine clinical immunoassays: Understand the causes and mitigate the risks. *Endocr Pract* 2017;23:989-98.
4. Li D, Radulescu A, Shrestha RT, Root M, Karger AB, Killeen AA, et al. Association of biotin ingestion with performance of hormone and nonhormone assays in healthy adults. *JAMA* 2017;318:1150-60.
5. Al-Salameh A, Becquemont L, Brailly-Tabard S, Aubourg P, Chanson P. A somewhat bizarre case of graves disease due to vitamin treatment. *Journal of the Endocrine Society* 2017;1:431-5.
6. Barbesino G. Misdiagnosis of graves' disease with apparent severe hyperthyroidism in a patient taking biotin megadoses. *Thyroid* 2016;26:860-3.
7. Bulow Pedersen I, Laurberg P. Biochemical hyperthyroidism in a newborn baby caused by assay interaction from biotin intake. *Eur Thyroid J* 2016;5:212-5.
8. De Roeck Y, Philipse E, Twickler TB, Van Gaal L. Misdiagnosis of graves' hyperthyroidism due to therapeutic biotin intervention. *Acta Clin Belg* 2018;73:372-6.
9. Evans N, Yates J, Tobin J, McGill J, Huynh T. Immunoassay interference secondary to therapeutic high-dose biotin: A paediatric case report. *J Paediatr Child Health* 2018;54:572-5.
10. Koehler VF, Mann U, Nassour A, Mann WA. Fake news? Biotin interference in thyroid immunoassays. *Clin Chim Acta* 2018;484:320-2.
11. Minkovsky A, Lee MN, Dowlatshahi M, Angell TE, Mahrokhian LS, Petrides AK, et al. High-dose biotin treatment for secondary progressive multiple sclerosis may interfere with thyroid assays. *AACE clinical case reports* 2016;2:e370-e3.
12. Stieglitz HM, Korpi-Steiner N, Katzman B, Mersereau JE, Styner M. Suspected testosterone- producing tumor in a patient taking biotin supplements. *J Endocr Soc* 2018;2:563-9.
13. Waghray A, Milas M, Nyalakonda K, Siperstein A. Falsely low parathyroid hormone secondary to biotin interference: A case series. *Endocr Pract* 2013;19:451-5.
14. FDA Safety Communication, November 28, 2017. <https://www.fda.gov/medical-devices/safety-communications/fda-warns-biotin-mayinterfere-lab-tests-fda-safety-communication>.
15. Favresse J, Burlacu MC, Maiter D, Gruson D. Interferences With Thyroid Function Immunoassays: Clinical Implications and Detection Algorithm. *Endocr Rev.* 2018 Oct 1;39(5):830-850.
16. Favresse J, Lardinois B, Nassogne MC, Preumont V, Maiter D, Gruson D. Anti-streptavidin antibodies mimicking heterophilic antibodies in thyroid function tests. *Clin Chem Lab Med.* 2018 Jun 27;56(7):e160-e163 *Clin Chem Lab Med.* 56(7): e160-e163.
17. Lam L, Bagg W, Smith G, Chiu WW, Middleditch MJ, Lim JC, Kyle CV. Apparent Hyperthyroidism Caused by Biotin-Like Interference from IgM Anti-Streptavidin Antibodies. *Thyroid.* 2018 Aug;28(8):1063-1067.
18. Harsch IA, Konturek PC, Böer K, Reinhöfer M. Implausible elevation of peripheral thyroid hormones during therapy with a protein supplement. *Clin Chem Lab Med.* 2017 Aug 28;55(9):e197-e198.
19. Peltier L, Massart C, Moineau MP, Delhostal A, Roudaut N. Anti-streptavidin interferences in Roche thyroid immunoassays: a case report. *Clin Chem Lab Med.* 2016 Jan;54(1):e11-4.
20. Rulander NJ, Cardamone D, Senior M, Snyder PJ, Master SR. Interference from anti-streptavidin antibody. *Arch Pathol Lab Med.* 2013 Aug;137(8):1141-6.

CONTACT

Veravas, Inc.
Research and Development
3510 Hopkins Place North
Oakdale, MN 55128
Phone: 1.888.466.4166
Email: info@veravas.com