

Resources

LEARNING AND EDUCATION PLATFORM

Animation

For more guidance on how to use Fetal Pillow please visit:
safeob.com/animation

PLAY

Online Training

To use the 15-minute online training tool visit:
training.safeob.com

REGISTER

Download the App



Scan the QR code using your mobile device

fetalpillow.beezer.com

1. Insertion



Fold the baseplate along the short axis ensuring balloon is on the inside.
Tube attachment must be at superior end (pointing upwards)
Insert into the vagina

2. Placement



Ensure balloon surface is in contact with the fetal head
Push device posteriorly, towards the coccyx
Placement is similar to a posterior Ventouse / vacuum cup

3. Legs Flat



Place patients legs flat on the table before inflation

4. Inflation



Use the syringe to inflate using 180ml of sterile saline
Close tap after filling to prevent fluid from leaking

5. Deliver baby

Perform Caesarean section using standard technique

6. Deflation & removal

Once baby has been delivered open the tap to drain the fluid
Remove the device at the end of the procedure