

Product information

REF 421



REF 421 **XL**

TRACOE experc Set *vario* XL

including **REF 451-P**

extra-long, spiral-reinforced, with adjustable neck flange, low-pressure cuff and minimally traumatic insertion system

The systems for *vario* XL consist of two parts: the white guiding catheter with silicone sleeve and the orange inserter. The silicone sleeve bridges the gap in diameter between the end of the tube and the inserter. The insertion system stabilizes and facilitates guiding the tube. The inserter fits the tube precisely and has a stop at the 15 mm connector. Once the tube is in the trachea, the colored inserter, the white guiding catheter with the silicone sleeve, and the Seldinger guide wire are removed together as a single unit. The silicone sleeve flaps back, and the unit can be removed.

The Benefits

Dilatation is performed using the Ciaglia one-step technique. The minimally traumatic inserter smooths over the gap in diameter between the end of the tube and the inserter. The stop prevents the inserter from being inserted too far into the tube, which reduces the risk of injury. The separately available P-tubes can also be used for reinsertion and/or for tube change – in combination with a Seldinger guide wire, if necessary.

experc Sets *vario*

TRACOE *vario* P-tubes REF 450-P, REF 451-P, REF 460-P, REF 461-P, REF 470-P and REF 471-P as well as the TRACOE experc Dilatation Set REF 520 are also available separately.

You can find details on the tube in the corresponding product information.

Accessories

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-CT
- ◆ TRACOE humid assist III REF 641
- ◆ TRACOE humid assist VI REF 642

For the Cuff Pressure Management

- ◆ TRACOE smart Cuff Manager REF 730-5
- ◆ TRACOE cuff pressure monitor (cpm) REF 720
- ◆ TRACOE cuff pressure monitor sensitive (cpms) REF 721