

LAB MONITORING SYSTEM

COMMON USE CASES + SENSORS

Freezers + refrigerators

- Temperature
- Electric power
- Door open/close
- Dry contact

Incubators

- Temperature
- Humidity
- CO2
- Door open/close
- Dry contact

Doors left open

- Freezers
- Refrigerators
- Incubators

Animal studies

- Light
- Temperature
- Humidity

Environmental parameters

- Light
- Temperature
- Humidity

Gas manifolds

- Dry contact
- System data port output

Also - ask us about pH Neutralization, Ovens, and more
