



Precision-Fit Owner's Manual

Cleaning Instructions

Precision-Fit seat cushions have two sewn removable incontinent lycra covers and Precision-Fit back cushions have one cover that can be machine washed, in cold water on the delicate cycle. Hang to dry.

If you have any questions regarding your cushion contact the vendor who provided it.

Mounting to Wheelchair

PRM Inc. does not attach the seating system to the wheelchair. Therefore, it is the responsibility of the vendor who provided your system to properly attach the seating to the wheelchair. Mounting hardware should be checked by the end user daily to ensure hardware is securely attached. Discontinue use if hardware has loosened and/or not properly attached. If you have any questions regarding your cushion contact the vendor who provided it.