























TRACOE Tracheostomy Tubes with Colour Coding

TRACOE *vario* tracheostomy tubes

TRACOE <i>vario</i>			TRACOE <i>vario</i> XL		
REF 450	spiral-reinforced, cuffed		REF 451 XL	spiral-reinforced, cuffed, extra-long	
REF 460	transparent, cuffed		REF 461 XL	transparent, cuffed, extra-long	
REF 470	transparent, cuffed, with subglottic suction line		REF 471 XL	transparent, cuffed, with subglottic suction line, extra-long	
REF 455	spiral-reinforced, without cuff				
REF 464	transparent, without cuff, fenestrated				

TRACOE *twist* tracheostomy tubes

TRACOE <i>twist</i>			TRACOE <i>twist</i> plus ++		
REF 301	cuffed, non-fenestrated		REF 311 ++	cuffed, non-fenestrated	
REF 302	cuffed, fenestrated		REF 312 ++	cuffed, double-fenestrated	
REF 303	without cuff, non-fenestrated		REF 313 ++	without cuff, non-fenestrated	
REF 304	without cuff, fenestrated		REF 314 ++	without cuff, double-fenestrated	
REF 306	cuffed, non-fenestrated, with subglottic suction line		REF 316 ++	cuffed, non-fenestrated, with subglottic suction line	
REF 888-306	cuffed, fenestrated, with subglottic suction line		REF 888-316 ++	cuffed, double-fenestrated, with subglottic suction line	
REF 305	Laryngectomy Tube, without cuff, non-fenestrated				
REF 309	cuffed, non-fenestrated, with air supply line for speaking				



















TRACOE Tracheostomy Tubes with Colour Coding

TRACOE *kids* tracheostomy tubes

TRACOE mini

REF 350	for neonates, without H ₂ O Cuff	 	REF 355	for children, without H ₂ O Cuff	 
----------------	--	---	----------------	--	---

TRACOE silcosoft®








REF 360	for neonates and infants, without H ₂ O Cuff, spiral-reinforced	 	REF 370	for children, without H ₂ O Cuff, spiral-reinforced	 
REF 361 PL	for neonates and infants, without H ₂ O Cuff, proximal longer, spiral-reinforced	 	REF 371 PL	for children, without H ₂ O Cuff, proximal longer, spiral-reinforced	 
REF 362	for neonates and infants, with H ₂ O Cuff, spiral-reinforced	 	REF 372	for children, with H ₂ O Cuff, spiral-reinforced	 
REF 363 PL	for neonates and infants, with H ₂ O Cuff, proximal longer, spiral-reinforced	 	REF 373 PL	for children, with H ₂ O Cuff, proximal longer, spiral-reinforced	 








TRACOE Tracheostomy Tubes without Colour Coding

TRACOE *twist* tracheostomy tubes

TRACOE *twist* inner cannulas





REF 401	inner cannulas, fenestrated, with speaking valve, pack of 2		REF 503	inner cannulas for laryngectomy tube, non-fenestrated, with 15 mm connector, pack of 3	
REF 501	inner cannulas, non-fenestrated, with 15 mm connector, pack of 3		REF 506	inner cannulas, fenestrated, with 15 mm connector, pack of 3	
REF 501-X	inner cannulas, non-fenestrated, with 15 mm connector, pack of 10		REF 506-X	inner cannulas, fenestrated, with 15 mm connector, pack of 10	
REF 501-XXX	inner cannulas, non-fenestrated, with 15 mm connector, pack of 30				

TRACOE *twist plus ++* inner cannulas

REF 521	inner cannulas, non-fenestrated, with 15 mm connector, pack of 3		REF 523	inner cannulas, double-fenestrated, with 15 mm connector, pack of 3	
REF 521-X	inner cannulas, non-fenestrated, with 15 mm connector, pack of 10		REF 523-X	inner cannulas, double-fenestrated, with 15 mm connector, pack of 10	
REF 521-XXX	inner cannulas, non-fenestrated, with 15 mm connector, pack of 30				

TRACOE *larynx* silicone short tubes


























TRACOE *larynx*

REF 580	Silicone Short Tube, with 22 mm housing, non-fenestrated, length 36 mm		REF 581	Silicone Short Tube, with 22 mm housing, fenestrated, length 36 mm	
REF 582	Silicone Short Tube, with 22 mm housing, non-fenestrated, length 53 mm		REF 583	Silicone Short Tube, with 22 mm housing, fenestrated, length 53 mm	



TRACOE Tracheostomy Tubes without Colour Coding

TRACOE *comfort* tracheostomy tubes

TRACOE <i>comfort</i>			TRACOE <i>comfort</i> XL / <i>comfort</i> plus ++		
REF 101	without cuff, non-fenestrated, without inner cannula		REF 201 XL	without cuff, non-fenestrated, extra-long, without inner cannula	
REF 102	without cuff, non-fenestrated, with inner cannula		REF 202 XL	without cuff, non-fenestrated, extra-long, with inner cannula	
REF 102-D	without cuff, non-fenestrated, with inner cannulas		REF 202-D XL	without cuff, non-fenestrated, extra-long, with inner cannulas	
REF 103	without cuff, fenestrated, with slide-on speaking valve type A, with 2 inner cannulas		REF 203 XL	without cuff, fenestrated, extra-long, with slide-on speaking valve type A, with 2 inner cannulas	
REF 103-A	without cuff, fenestrated, with slide-on speaking valve type A, with 3 inner cannulas		REF 203-A XL	without cuff, fenestrated, extra-long, with slide-on speaking valve type A, with 3 inner cannulas	
REF 104	without cuff, fenestrated, with swivel speaking valve type B, with 2 inner cannulas		REF 204 XL	without cuff, fenestrated, extra-long, with swivel speaking valve type B, with 2 inner cannulas	
REF 104-A	without cuff, fenestrated, with swivel speaking valve type B, with 3 inner cannulas		REF 204-A XL	without cuff, fenestrated, extra-long, with swivel speaking valve type B, with 3 inner cannulas	
REF 105	without cuff, non-fenestrated, with 15 mm connector		REF 205 XL	without cuff, non-fenestrated, extra-long, with 15 mm connector	
REF 105-D	without cuff, non-fenestrated, with 2 x 15 mm connector				
REF 115	without cuff, fenestrated, with 15 mm connector		REF 215 XL	without cuff, fenestrated, extra-long, with 15 mm connector	
REF 106-D	Laryngectomy Tube, without cuff, non-fenestrated, with 2 inner cannulas				
REF 107	without cuff, fenestrated, with 15 mm connector, with 22 mm housing		REF 157 ++	without cuff, double-fenestrated, with maximised length, with 15 mm connector, with 22 mm housing	
REF 108	without cuff, fenestrated, with swivel speaking valve type B, with 22 mm housing		REF 158 ++	without cuff, double-fenestrated, with maximised length, with swivel speaking valve type B, with 22 mm housing	
REF 114	without cuff, fenestrated, with swivel speaking valve type B, with 15 mm connector		REF 164 ++	without cuff, double-fenestrated, with maximised length, with swivel speaking valve type B, with 15 mm connector	