

Product information

REF 580 and 582



TRACOE *larynx* Silicone Short Tube with 22 mm housing

REF 580


Length: ◆ 36 mm
Including: ◆ neck strap

REF 582

Length: ◆ 53 mm
Including: ◆ neck strap

Silicone short tubes are ideal for the use in stable tracheostomas to secure the respiratory tract once no tracheostomy tube is needed anymore.

Technical Data



REF 580, REF 581, REF 582, REF 583

Size	ID mm	OD mm	REF 580 REF 581		REF 582 REF 583	
			TL mm	OB mm	TL mm	OB mm
08	9.5	12.0	36	43	53	62
09	10.5	13.5	36	43	53	62
10	12.0	15.0	36	43	53	62
12	13.5	17.0	36	43	53	62


ID: inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;
OB: length along outer bend; **L:** length

To order, specify **REF + size**, e.g. **REF 601-08**

Product information

REF 580 and 582

Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Accessories

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S

TRACOE larynx Silicone Short Tube is available in lengths of 36 mm with order number REF 580 and by order number REF 582 of 53 mm.