

Product information

REF 581 and 583



TRACOE *larynx* Silicone Short Tube
with 22 mm housing, fenestrated

REF 581



Length: ◆ 36 mm
Including: ◆ neck strap


REF 583



Length: ◆ 53 mm
Including: ◆ neck strap

For tracheostomised patients as well as laryngectomised patients with a voice prosthesis.

Technical Data



REF 580, REF 581, REF 582, REF 583

Size	ID mm	OD mm	REF 580 REF 581		REF 582 REF 583	
			TL mm	OB mm	TL mm	OB mm
08	9.5	12.0	36	43	53	62
09	10.5	13.5	36	43	53	62
10	12.0	15.0	36	43	53	62
12	13.5	17.0	36	43	53	62


ID: inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;
OB: length along outer bend; **L:** length

To order, specify **REF + size**, e.g. **REF 601-08**

Product information

REF 581 and 583

Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Accessories

Speaking Valves and Stopper

- ◆ TRACOE phon assist I REF 650-S, REF 650-SO
- ◆ TRACOE phon assist II REF 655-S
- ◆ TRACOE modular stopper REF 621

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S

TRACOE larynx Silicone Short Tube with fenestration is available in lengths of 36 mm with order number REF 581 and by order number REF 583 of 53 mm.