

Product information

REF 612 and 613



TRACOE larynx stoma button
made of silicone

REF 612


Length: ♦ 22mm
♦ in sterile package
Sizes: ♦ 06 – 12

REF 613

Length: ♦ 30mm
♦ in sterile package
Sizes: ♦ 06 – 12

Technical Data

REF 601, REF 602, REF 603, REF 612, REF 613



Size	ID mm	OD mm	REF 601	REF 602 REF 612	REF 603 REF 613
			short	standard	long
			L mm	L mm	L mm
06	6	9	16	22	30
07	7	10	16	22	30
08	8	11	16	22	30
09	9	12	16	22	30
10	10	13	16	22	30
11	11	14	16	22	30
12	12	15	16	22	30


ID: inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;
OB: length along outer bend; **L:** length

To order, specify **REF + size**, e. g. **REF 601-08**

Product information

REF 612 and 613

Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Accessories

Speaking Valves and Stopper

- ◆ TRACOE phon assist I REF 650-S, REF 650-SO
- ◆ TRACOE phon assist II REF 655-S
- ◆ TRACOE modular stopper REF 621

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S