

# Product information

## REF 601



### REF 601

**TRACOE grid button**  
short, made of silicone


Including: ♦ neck strap  
♦ in sterile package

Sizes: ♦ 06 – 12

TRACOE grid button with eyes for attaching a neck strap and with a removable grid to prevent inhalation of foreign objects.

### Technical Data

REF 601, REF 602, REF 603, REF 612, REF 613



Size	ID mm	OD mm	REF 601	REF 602 REF 612	REF 603 REF 613
			short	standard	long
			L mm	L mm	L mm
06	6	9	16	22	30
07	7	10	16	22	30
08	8	11	16	22	30
09	9	12	16	22	30
10	10	13	16	22	30
11	11	14	16	22	30
12	12	15	16	22	30


**ID:** inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;  
**OB:** length along outer bend; **L:** length

To order, specify **REF + size**, e. g. **REF 601-08**

# Product information

## REF 601

### Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

### Accessories

#### Speaking Valves and Stopper

- ◆ TRACOE phon assist I REF 650-S, REF 650-SO
- ◆ TRACOE phon assist II REF 655-S
- ◆ TRACOE modular stopper REF 621

#### Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S