

Product information

REF 602

REF 602

TRACOE grid button
standard, made of silicone

Including: ♦ neck strap
♦ in sterile package


Sizes: ♦ 06 – 12



TRACOE grid button with eyes for attaching a neck strap and with a removable grid to prevent inhalation of foreign objects.

Technical Data

REF 601, REF 602, REF 603, REF 612, REF 613



Size	ID mm	OD mm	REF 601	REF 602 REF 612	REF 603 REF 613
			short	standard	long
			L mm	L mm	L mm
06	6	9	16	22	30
07	7	10	16	22	30
08	8	11	16	22	30
09	9	12	16	22	30
10	10	13	16	22	30
11	11	14	16	22	30
12	12	15	16	22	30


ID: inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;
OB: length along outer bend; **L:** length

To order, specify **REF + size**, e. g. **REF 601-08**

Product information

REF 602

Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Accessories

Speaking Valves and Stopper

- ◆ TRACOE phon assist I REF 650-S, REF 650-SO
- ◆ TRACOE phon assist II REF 655-S
- ◆ TRACOE modular stopper REF 621

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S