

Product information

REF 603



REF 603

TRACOE grid button
long, made of silicone


Including: ♦ neck strap
♦ in sterile package

Sizes: ♦ 06 – 12

TRACOE grid button with eyes for attaching a neck strap and with a removable grid to prevent inhalation of foreign objects.

Technical Data

REF 601, REF 602, REF 603, REF 612, REF 613



| Size | ID mm | OD mm | REF 601 | REF 602 REF 612 | REF 603 REF 613 |
|------|----------|----------|---------|--------------------|--------------------|
| | | | short | standard | long |
| | | | L mm | L mm | L mm |
| 06 | 6 | 9 | 16 | 22 | 30 |
| 07 | 7 | 10 | 16 | 22 | 30 |
| 08 | 8 | 11 | 16 | 22 | 30 |
| 09 | 9 | 12 | 16 | 22 | 30 |
| 10 | 10 | 13 | 16 | 22 | 30 |
| 11 | 11 | 14 | 16 | 22 | 30 |
| 12 | 12 | 15 | 16 | 22 | 30 |

ID: inside diameter; **OD:** outside diameter of tube; **TL:** length along centre line;
OB: length along outer bend; **L:** length

To order, specify **REF + size**, e. g. **REF 601-08**

Product information

REF 603

Scope of Supply of TRACOE *larynx* Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Accessories

Speaking Valves and Stopper

- ◆ TRACOE phon assist I REF 650-S, REF 650-SO
- ◆ TRACOE phon assist II REF 655-S
- ◆ TRACOE modular stopper REF 621

Heat/Moisture Exchangers, HME

- ◆ TRACOE humid assist I REF 640-S