

Product information

REF 607-45 and 607-70



TRACOE button plus

Silicone decannulation button made of silicone

REF 607-45


- Length: ◆ 45 mm
- ID: ◆ 8 mm
- OD: ◆ 11 mm
- Including: ◆ neck strap
- ◆ in sterile package

REF 607-70

- Length: ◆ 70 mm
- ID: ◆ 8 mm
- OD: ◆ 11 mm
- Including: ◆ neck strap
- ◆ in sterile package

For the therapy of tracheostomised patients after long-term ventilation – safe airway after decannulation. The sealing is very airtight and dislocation save. The patient is able to cough up the stimulus secretion very fast. With the adjustable neck flange the TRACOE button plus can be positioned optimal to the soft neck tissue.

Scope of Supply of TRACOE larynx Products



All TRACOE larynx silicone short tubes, grid buttons and button plus products are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

Caution, not suitable for:

- ◆ laryngectomised patients
- ◆ patients with unstable breathing
- ◆ patients, who primarily breathe through the tracheostoma

TRACOE button plus is available in lengths of 45 mm with order number REF 604-45 and by order number REF 607-70 of 70 mm.