

# Product information

## REF 106-D




### REF 106-D

#### TRACOE *comfort* Laryngectomy Tube

Including: ♦ 2 inner cannulas with low profile  
♦ neck strap

For patients whose larynx has been removed.

### Technical Data



#### REF 106-D

Size	ID-OC mm	OD-OC bc mm	OD-OC bf mm	ID-IC mm	TL mm	OB mm	BA °
07	7.0	8.8	10.5	5.8	47	50	130°
08	8.0	10.0	11.3	6.8	50	53	130°
09	9.0	11.4	12.6	7.6	52	57	130°
10	10.0	12.4	13.7	8.7	55	60	130°
11	11.0	13.1	15.0	9.7	58	63	130°
12	12.0	14.2	16.2	10.8	60	66	130°


**ID-OC:** inside diameter (clear width) at bottom of outer cannula; **OD-OC bc:** outside diameter at bottom of outer cannula;  
**OD-OC bf:** outside diameter of outer cannula behind the flange; **ID-IC:** inside diameter at bottom of inner cannula;  
**TL:** length along centre line from start of neck flange to bottom of tube; **OB:** length along outer bend (greatest length);  
**BA:** bending angle

To order, specify **REF + size**, e.g. **REF 103-08**

# Product information

## REF 106-D

### Scope of Supply of TRACOE *comfort* Tubes



All TRACOE *comfort* tubes are supplied with a wide adjustable neck strap and a product info card with two detachable labels.

### Cleaning products

- ◆ TRACOE care Liquid Paraffin REF 902  
lets inner cannulas be easily inserted and removed, assures best possible fit and good quality of speech
- ◆ TRACOE tube clean Cleaning Kit REF 930-A/-B  
with cleaning powder REF 932
- ◆ TRACOE tube clean Cleaning Kit REF 931-A/-B  
with concentrated cleaning liquid REF 933

### Product Variants

Product variants can be made in order to optimise the functionality and wearing comfort and also for patients with complex anatomical and physiological conditions. The relevant ordering forms are available at [www.tracoe.com/en/downloads](http://www.tracoe.com/en/downloads).