

Product information

REF 850


REF 850

TRACOE Pleural Drainage Set using the Seldinger Technique

The seldinger pleural catheter is used for draining pleural effusion and to relieve a pneumothorax.

Packaging unit:
Including:

- ◆ 5 units
- ① drainage catheter with inner stylet, radiopaque, kink-resistant due to plastic spiral reinforcement, MR conditional, luer lock connector
- ② scalpel (11), single use
- ③ syringe (10ml)
- ④ tuohy puncture needle (1.2 x 1.5 x 90mm, Gauge 17)
- ⑤ coated seldinger guide wire with Nitinol core
- ⑥ dilator with depth limiter
- ⑦ fixation plaster
- ⑧ 3-way-valve
- ⑨ adapter from luer lock to funnel connector
- ⑩ secretion bag (1.5l) with luer lock connection



The set's features include:

- ◆ kink-resistant, thin drainage catheter (12 Fr or 14 Fr) with embedded plastic spiral reinforcement
- ◆ extra-soft catheter tip with three staggered openings
- ◆ tuohy puncture needle
- ◆ seldinger guide wire, kink-resistant due to Nitinol core
- ◆ convenient thumb inserter
- ◆ maximum in situ usage time of 28 days
- ◆ radiopaque and MR conditional