


Product information

REF 520



REF 520

TRACOE experc Dilatation Set for percutaneous tracheostomy

- 1 scalpel
- 2 syringe
- 3 puncture needle (14G) with teflon catheter
- 4 seldinger guide wire made of kink resistant Nitinol with inserter
- 5 short dilator (14Ch/Fr)
- 6 guiding catheter with safety stop
- 7 TRACOE experc dilator with hydrophilic coating (becomes slidable after moistening)
- 8 4 compresses

◆ in sterile package

The TRACOE experc Dilatation Set (REF 520) can be used alone or together with a P-tube most appropriate for the patient.