



- Contents of 2 Covid-19 Self Tests:**
- 2 On/Go COVID-19 test cassettes
 - 2 Nasal swabs
 - 2 Reagent vials
 - 2 Vial caps

Store between 33.8-86°F (1-30°C) until use

If you have symptoms of COVID-19, you can use a single test. If you do not have symptoms of COVID-19, you will need at least two tests per person. You may need to purchase additional tests to perform serial (repeat) testing. This test is more likely to give you a false negative result when you have COVID-19 than a lab-based molecular test.

This product has not been FDA cleared or approved, but has been authorized by FDA under an EUA;

This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564 (b) (1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.



Distributed by:
Intrivo
2021 Santa Monica Blvd, #11
Santa Monica, CA 90404, USA
888-965-0301
www.intrivo.com



Manufactured by:
Access Bio, Inc.
65 Clyde Road, Suite A,
Somerset, NJ 08873, USA
732-873-4040
www.accessbio.net



on/go

For Use Under Emergency Use Authorization

on/go


**COVID-19
Antigen Self-Test**

**RESULTS
IN 10 MIN**

**SMART
APP**

**2
SELF
TESTS**

iOS13 or newer
for Apple iPhone
Android10 or newer
for Android Phone



Scan to download
the On/Go™ mobile
app to start test.
www.letsongo.com

For use under
Emergency Use
Authorization
(EUA) only. For in
vitro diagnostic
use only. For use
with lower nasal
swab specimens.



- Contents of 1 Covid-19 Self Test:**
- 1 On/Go COVID-19 test cassette
 - 1 Nasal swab
 - 1 Reagent vial
 - 1 Vial cap

Store between 33.8-86°F (1-30°C) until use

If you have symptoms of COVID-19, you can use a single test. If you do not have symptoms of COVID-19, you will need at least two tests per person. You may need to purchase additional tests to perform serial (repeat) testing. This test is more likely to give you a false negative result when you have COVID-19 than a lab-based molecular test.

This product has not been FDA cleared or approved, but has been authorized by FDA under an EUA;

This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and,

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564 (b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Intrivo

Distributed by:
Intrivo
2021 Santa Monica Blvd, #11
Santa Monica, CA 90404, USA
888-965-0301
www.intrivo.com


AccessBio

Manufactured by:
Access Bio, Inc.
65 Clyde Road, Suite A,
Somerset, NJ 08873, USA
732-873-4040
www.accessbio.net

UDI Placement



For Use Under
Emergency Use
Authorization



REF: RPM-00179

on/go

on/go

on/go™

on/go

**COVID-19
Antigen Self-Test**

RESULTS
IN 10 MIN

iOS13 or newer
for Apple iPhone
Android10 or newer
for Android Phone

SMART
APP

Download on the
App Store
GET IT ON
Google Play

1 SELF
TEST



Scan to download
the On/Go™ mobile
app to start test.
www.letsongo.com

For use under
Emergency Use
Authorization
(EUA) only. For in
vitro diagnostic
use only. For use
with lower nasal
swab specimens.