

Dayspring®

MEASUREMENT FORM (1/2)

DATE

Patient Information

PATIENT NAME

AFFECTED LIMB
(CIRCLE ONE)

PATIENT DATE OF BIRTH

LEFT | RIGHT | BOTH

Person Performing Measurements

NAME

TITLE

FACILITY / LOCATION

Shipping

NAME

COMPANY

STREET ADDRESS

CITY

STATE / PROVINCE

ZIP / POSTAL CODE PHONE

EMAIL

(FOR SHIPPING CONFIRMATION)

General Measurements


Referring to the diagram provided at right, complete the corresponding fields below.
All measurements should be taken in centimeters (cm) to the nearest half unit.

MEASUREMENT (A) DISTANCE FROM WRIST TO ELBOW: cm

MEASUREMENT (B) DISTANCE FROM WRIST TO AXILLA: cm

MEASUREMENT (C) DISTANCE FROM WRIST TO KNUCKLES: cm

MEASUREMENT (D) CIRCUMFERENCE OF KNUCKLES: cm



Arm Garment	Arm Measurement B (wrist to axilla)	Arm Small F-0010-01	Less than 44cm	Max circumference at Axilla 52cm
		Arm Med F-0010-02	Between 44cm & 48cm	
		Arm Long F-0010-03	Greater than 48cm	

Proceed to Page 2.

Dayspring®

MEASUREMENT FORM (2/2)

Circumference Measurements


The size of the garment depends on the length of the patient's arm (B), measured previously on page one. Mark and complete the corresponding section below.

Short:
(B) is below 44cm

DISTANCE FROM WRIST 0 4 11 18 25 31 38
cm cm cm cm cm cm cm

CIRCUMFERENCE

cm	cm	cm	cm	cm	cm	cm	cm
----	----	----	----	----	----	----	----




SECTION 4S

Medium:
(B) is between
44cm and 48cm

DISTANCE FROM WRIST 0 4 12 20 27 35 43
cm cm cm cm cm cm cm

CIRCUMFERENCE

cm	cm	cm	cm	cm	cm	cm	cm
----	----	----	----	----	----	----	----




SECTION 4M

Long:
(B) is above
48cm

DISTANCE FROM WRIST 0 4 12 19 26 33 41 48
cm cm cm cm cm cm cm

CIRCUMFERENCE

cm	cm	cm	cm	cm	cm	cm	cm
----	----	----	----	----	----	----	----



SECTION 4L