

Centricity Universal Viewer

Zero Footprint Client

Version 6.0, SP4

Features

- No software download to the desktop
- Zero administrative rights required on the device
- Uses HTML5 standards
- Desktop version for diagnostic purposes
- Mobile devices for review only mode (not for primary diagnosis)
- Support for dual monitor configuration (requires Native Client Websocket installation)

Supported Browsers

ZFP desktop

- Internet Explorer 10 and 11 (32 bit and 64 bit)
- Internet Explorer 8 and 9 (32 bit) with the Legacy Browser Component
- Microsoft Edge (64 bit) version 10525
- Mozilla Firefox version 44
- Google Chrome version 48 on Windows
- Google Chrome version 48 on Macintosh
- Apple Safari version 8 on Macintosh
- ✓ **Note:** For the browser instance that is running ZFP DICOM Viewer, there should be 1024 MB minimum memory when using Internet Explorer 9, and 512 MB minimum memory available for other supported browsers.
- ✓ **Note:** For the browser instance that is running XDS-Based Patient Timeline, there should be 2 GB minimum memory available when using all supported browsers.

ZFP mobile devices

- Apple Safari iOS 8 and iOS 9 on iPad with Retina Display
- Google Chrome version 48 with iOS 9 on iPad
- Google Chrome version 48 on Android KitKat v4.4.2
- ✓ **Note:** Google Chrome version 45 with iOS 9 cannot play .MOV or .MPEG files. Download the files to an external application

Gaining Access

- Visit zfp.amimaging.com and enter your account credentials at the login page
- If you do not have an account, you can submit a request by one of the following methods:
 - ✓ Visit www.amimaging.com, click the Physician Login link in the top right corner, and choose "System Access Request Form." Download the form and follow the instructions to return it.
 - ✓ Request a form by emailing RPACS@amimaging.com or calling 402-484-6677