

# Philips HeartStart Automated External Defibrillator

**Simple to operate  
with icon-driven  
and audible  
guided  
instructions come  
together in this  
intuitive AED.**



## **Guided audio instruction talks you through your training**

The HeartStart OnSite provides on-demand audio coaching. Calm, clear voice instructions are tailored to the responder's actions, providing AED and CPR guidance during the resuscitation of a sudden cardiac arrest victim.

## **Ready when needed**

Your HeartStart OnSite goes through a 120-point test before it leaves the factory. On the job, it conducts self tests every single day, not just weekly. It performs over 80 different tests in all. Even the pads are tested for readiness.

## **Quick Shock feature**

Minimize chest compression interruptions CPR is even more vital to survival than previously realized. Rapidly delivering a shock after chest compressions is critical. The HeartStart OnSite's Quick Shock feature reduces the time between hands-off and shock delivery.

## **Integrated SMART pads**

SMART pads sense application to patient skin and the pads are "smart" because they can sense that they have been removed from the cartridge and applied to the patient. This helps the HeartStart OnSite know at what point the user is in the response process.