

Arrow

# Intra-Aortic 50 cc Balloon Catheters

## With New Sizing Guidelines

### Meeting your patients' counterpulsation needs

#### Sizing Guidelines\*

Arrow IAB Catheters	30 cc	40 cc	50 cc
Height (feet/inches)	4' 10" - 5' 4"	5' 4" - 6'	> 6'
Height (centimeters)	147 - 162 cm	162 - 182 cm	> 182 cm
Body Surface Area	< 1.8 m <sup>2</sup>	> 1.8 m <sup>2</sup>	> 1.8 m <sup>2</sup>

\*Refer to IFU

#### 50 cc FiberOptix IAB Catheters

- FiberOptix Sensor delivers high fidelity AP signal and enables proprietary IntraBeat Timing when connected to AutoCAT2WAVE and AC3 Optimus IABP Systems
- 0.027" thin-walled alloy central lumen with 1.5" flex tip design
- Coil-reinforced central lumen to help reduce potential for catheter whip
- Robust polycarbonate bifurcate with large bore connector ensures high rate helium exchange up to 180 bpm

#### 50 cc UltraFlex IAB Catheters

- The only coil reinforced 50 cc IAB catheters currently available
- Designed with features:
  - Navigating tortuous anatomy
  - Flexibility
  - Kink resistance to promote optimal gas exchange
- 0.027" thin-walled alloy central lumen with 1.5" flex tip design
- Coil reinforced catheter to help reduce potential for catheter whip
- Robust polycarbonate bifurcate with large bore connector ensures high rate helium exchange up to 180 bpm

#### 50 cc RediGuard IAB Catheters

- Flexible, soft catheter with stainless-steel braid central lumen
- Centimeter markings on catheter body to facilitate placement